



# BASIC

## Defibrillator



Turkish First Aid Association supports "Life-Point"

CE  
1984

# BASIC

## Defibrillator

### Technical Specifications

Waveform	Monophasic Waveform
Operating Modes	Manual / Synchronous
Energy Selection (Joule)	5-10-30-50-100-150-200-300-360
Working Conditions	0 / 40 °C and % 30-95 Rh
Storing Conditions	-20 / +55 °C
Dimensions	40cm x 15cm x 33cm
Weight	7 Kg (included battery and paddles)
Main Supply Voltage	AC 220V 50-60 Hz ( $\pm$ % 20)
Pediatric Paddles	Yes
Trolley	Optional
Ambulance Bracket (EN1789)	Optional

\* METsis Medikal has right to amend the values written above any time.



**METsis**  
ADVANCED TECHNOLOGY SYSTEMS



**METUTECH**  
ODTÜ-TEKNOKENT

Ar-Ge | R&D  
ODTÜ Ostim Teknokent İdari Binası  
Uzay Çağı Caddesi 78. Sokak No: 6  
2.B Kat D: 1-2 Ostim-Ankara/Turkey  
Tel: +90.312 386 21 90  
Faks: +90.312 386 21 80

info@metsismedikal.com

