




PRO Biphasic

Defibrillator



Turkish First Aid Association supports "Life-Point"

CE
1984

PRO Biphasic

Defibrillator

Technical Specifications

Waveform	Biphasic Waveform
Operating Modes	Manual / Synchronous
Energy Selection (Joule)	2-5-7-10-20-30-70-100-150-200-270
ECG Derivations	DI DII DIII AvR AvL AvF V1...V6
ECG Filters	Monitor I Monitor II Strong Surgery Diagnostic
Heart Rate Limits	0 - 300 bpm
Monitor Type	5,7 inch LCD
Monitor Speed	25mm/sec - 50mm/sec
Printer Type / Speed	Thermal Printer / 25 - 50 mm/sec
Printer Paper Size	50mm x 30m
Maximum Discharging (270 J)	100 (with full battery)
ECG Monitoring Time	4 hours (with full battery)
Working Conditions	0 / 40 °C and % 30-95 Rh
Storing Conditions	-20 / +55 °C
Dimensions	40cm x 15cm x 33cm
Weight	8,5 Kg (included battery and paddles)
Main Supply Voltage	AC 220V 50-60 Hz (± % 20)
Pediatric Paddles	Yes
SpO ₂	Optional
NIBP	Optional
Trolley	Optional
Ambulance Bracket (EN1789)	Optional

* METsis Medikal has right to amend the values written above any time.



METsis
ADVANCED TECHNOLOGY SYSTEMS



METUTECH
ODTÜ-TEKNOKENT

Ar-Ge | R&D

ODTÜ Ostim Teknokent İdari Binası

Uzay Çağı Caddesi 78. Sokak No: 6

2.B Kat D: 1-2 Ostim-Ankara/Turkey

Tel: +90.312 386 21 90

Faks: +90.312 386 21 80

info@metsismedikal.com

