



## Plasma Sterilizer

Sterilmed Medical Plasma Sterilizer is an Eco-Friendly, Low Temperature Sterilizer used in sterilizing thermo-sensitive hospital equipment and materials with the use of H<sub>2</sub>O<sub>2</sub> (Hydrogen peroxide).

### General Features

High Quality Stainless Steel (316 or 304 AISI),  
Used in hospitals, laboratories, prosthetics manufacturers, dental and veterinary,  
Easy to use Plug and Play device,  
Rapid Warmup & Drying System,  
Safe to use, yields no residue (only water and oxygen,)  
Fast running cycle to increase turnover rate of equipment and materials.

### Safe and Sound

With security measures that protect the operator and the device.  
Emergency stop button,  
Sensors ensure that doors are locked,  
Doors locks during procedure, Unable to open both doors at once in, Septic-Aseptic models.

### Control and Observe

Every step on a handy 7" Touchscreen and PLC,  
Full-automatic programs, adjustable for specific calibrations,  
User-Friendly, easy to use UI,  
Numeric and graphic display,  
Integrated 40-character thermal printer,

Visual, written and audial warning systems,  
Data storage of up 1000 cycles,  
USB port for transferring saved data to other media,  
Communications via RS 232 and RS485 Ethernet Module.



[Compatible Integrated Baskets](#)



[Compatible Hydrogen Peroxide Cartridge System](#)



### Versatility

Of our devices makes room for a variety of additional features and accessories such as:



Computer Connection Software



Gsm Connection Module



10" Color TFT LCD Touchscreen



Rack and Baskets



### For a Green Future

As Sterilmed Medical, we want to leave a green future for next generations and we manufacture our products with that sensibility; a zero waste mentality. All materials used in manufacturing our devices are recyclable or eco-friendly and we constantly work on increasing the efficiency of devices for less water and energy consumption.



[www.sterilmed.com.tr](http://www.sterilmed.com.tr)



### Quality Without Compromise

Is our motto and we manufacture our devices assuring their quality in compliance

Medical Devices Regulations:  
MDR 2017/745/EU

Medical Devices Class:  
Class IIa, acc. To EC MDR 2017/745/EU Annex IX

Low Voltage Directives:  
EN 60601-1

Electromagnetic Compatibility Directives:  
EN 60601-1-2

Sterilization of health care products:  
ISO 14937

Quality Management System Requirements:  
ISO 9001

Medical Devices Regulatory Requirements:  
ISO 13485

Environment Management System:  
ISO 14001

Model	Chamber (WxHxD) - mm			Litres	Overall (wxhxd) - mm		
PL02	400	400	650	100	750/680	1760/1640	820/740
PL03	400	400	800	125	750/680	1760/1640	980/870
PL04	400	400	950	150	750/680	1760/1640	1130/1020

\* Single / Double Automatic Vertical Sliding Door System and Manuel Hinged Door System options are available.

As the main manufacturer of STERILMED Products, besides Standard STERILMED Model products and designs, visual changes can be made on a project-specific basis.



www.sterilmed.com.tr



### Quality Certificates



Başkent O.S.B. | +90 312 375 81 00  
18. Cadde No.43 Malıköy, | info@sterilmed.com.tr  
Sincan, Ankara/Türkiye | www.sterilmed.com.tr



# Plasma Sterilizer and Equipments



www.sterilmed.com.tr