



Neon

Proximal Femoral Nail

- Nail length: Neon Short: 220 mm & Neon Long: 300 mm to 420 mm in 20 mm increment
- Proximal Diameter: 15.5 mm
- Distal Diameter: 10, 11 mm
- Lag Screw Angle: 120°, 125° and 130°




Ntan

Proximal Femoral Nail

- Nail length: 180, 200, 220 mm
- Proximal Diameter: 16.3 mm
- Distal Diameter: 10, 11.5, 13 mm
- Lag Screw Angle: 130°

Navy

A/R Femoral Nail



- Nail length: 280 mm to 460 mm in 20 mm increment
- Proximal Diameter: 13 mm
- Distal Diameter: 10, 11, 12, 13 mm
- Compression Range: 10 mm



Navy - S

Universal A/R Femoral Nail

Claw + Screw
One Instrument Set
Femur + Tibia

Nite

Tibial Nail



- Nail length: 255 mm to 375 mm in 15 mm increment
- Proximal Diameter: 11 mm
- Distal Diameter: 9, 10, 11, 12 mm
- Compression Range: 7 mm



Nite - S

Universal Tibial Nail



Nice

Canulated Compression Screw

- Length: 45 to 120 mm in 5 mm increment
- Shaft Diameter: 6.2 mm
- Distal Thread Diameter: 7.5 mm
- Distal Thread Length: 16 mm
- Compression Range: 4 - 8 mm



Nuna

Humeral Nail


- Nail length: 165 mm to 315 mm in 10 mm increment
- Proximal Diameter: 10 mm
- Distal Diameter: 8, 9, 10 mm

Lag Screw Claw Technology



One Instruments Set - Four IM. Nails

Femur Claw Antegrade
Femur Claw Retrograde
Femur Screw Antegrade
Femur Screw Retrograde



Tibia Claw
Tibia Screw

Six Surgical Options



- ✓ Increased stability
- ✓ Shorter operative time
- ✓ Lower radiation exposure
- ✓ Fewer incisions
- ✓ No free-hand locking

Ntan

Proximal Femoral Nail

- Nail length: 180, 200, 220 mm
- Proximal Diameter: 16.3 mm
- Distal Diameter: 10, 11.5, 13 mm
- Lag Screw Angle: 130°
- Lag Screw length: 70 - 120 mm with 5mm



Navy-X

Femoral Nail

- 3 options: Recon, Transverse, Oblique
- Nail length: 320 mm to 420 mm in 20 mm increment
- Distal Diameter: 10, 11, 12, 13 mm
- Recon Screws 6.4mm, 65mm to 115mm with 5mm increment


Nite-X Expert

Expert Tibial Nail


- Nail length: 285 mm to 360 mm in 15 mm increment
- Distal Diameter: 9, 10, 11 mm





Flexible Reamer Set 1



Flexible Reamer Set 2





Neon PFN Set



Ntan PFN Set



Femoral Nail Set Type II



Tibial Nail Set Type II

