



GEOTEK[®]
Healthcare Products

ABOUT US

Geotek that was founded in Ankara in 2003 produces single use medical consumables. As a result of increasingly continuing R&D and Product Development activities; it renews the production technology, increases range of products in its portfolio day by day and presents to service of end users.

Geotek has assessed dynamics within time elapsed from production of materials until arrival to end user analytically by working in continuous contact with health professionals and becomes a successful and reliable company by keeping product and company service quality at the top level.

We are in a responsibility that the pride of being a successful and reliable company should have continuity also in the future like the past. Therefore, we shall keep this continuity by proceeding to do our business carefully and grow together with you by protecting our success, confidence, care and stability values.

GEOTEK®
Healthcare Products


CONTENTS

BIOPSY NEEDLES

Katracore® Automatic Biopsy Needle	4
Vescut® Automatic Biopsy Needle	5
Estacore® Automatic Biopsy Needle	6
Alfa Gun® Reusable Biopsy Device	7
Maxicore® Reusable Biopsy Gun	8
Maxicore® Biopsy Gun Needle	9
Maxicore-M® Biopsy Gun Needle	10
Maxicore-U® Biopsy Gun Needle	11
Semi-Automatic Biopsy Needle NEW DESIGN	12
Semi-Automatic Biopsy Needle	13
Chiba Biopsy Needle	14
Coaxial Needle	15
Breast Localization Needle	16
Percutaneous Access Needle	17
Tumor Marking Needle w/PEEK Marker	18
Tumor Marking Needle w/Gold Marker	19
Bone Marrow Aspiration Needle	20
Bone Marrow Biopsy Needle	21
Disposable Needle Guides	22
Trus Biopsy Kit	23
Oocyte Aspiration Needle	24

UROLOGY CATHETERS

Long Durable DJ Ureteral Stent Set	26
Double J Ureteral Stent Set	27
Ureteral Catheter	28
Single J Stent Set	29
Re-entry Malecot Catheter	30
Percutaneous Nephrostomy Set	31
Direct Puncture Nephrostomy Set	32
Ureteral Balloon Dilator	33
Suprapubic Drainage Catheter Set	34
Amplatz Renal Dilator Set	35
Stone Basket	36
Stone Grasper	37
Nephrostomy Balloon Dilator	38
Ureteral Access Sheath	39




BIOPSY NEEDLES

KATRACORE® AUTOMATIC BIOPSY NEEDLE

- ↳ Easy and one hand operated with two stage arming
- ↳ Echogenic and atraumatic needle tip design
- ↳ 22mm advancement with 20mm notch
- ↳ Ergonomic and light weight design
- ↳ Depth marked canula
- ↳ Coaxial needles are optionally available



Ref. No	Diameter (G)	Length (cm)	Coaxial Needle Length (cm)	With Coaxial Ref. No
GKR1410	14	10	5,5	GKR1410-X
GKR1415	14	15	10,5	GKR1415-X
GKR1420	14	20	15,5	GKR1420-X
GKR1610	16	10	5,5	GKR1610-X
GKR1615	16	15	10,5	GKR1615-X
GKR1620	16	20	15,5	GKR1620-X
GKR1810	18	10	5,5	GKR1810-X
GKR1815	18	15	10,5	GKR1815-X
GKR1820	18	20	15,5	GKR1820-X
GKR1825	18	25		
GKR1830	18	30		
GKR2010	20	10	5,5	GKR2010-X
GKR2015	20	15	10,5	GKR2015-X
GKR2020	20	20	15,5	GKR2020-X



VES CUT® AUTOMATIC BIOPSY NEEDLE

- ↳ Choice of two firing buttons
- ↳ Easy and one hand operated with two stage arming
- ↳ Automatically activated safety mechanism
- ↳ Echogenic and atraumatic needle tip design
- ↳ Ergonomic and light weight design
- ↳ Depth marked canula
- ↳ Coaxial needles are optionally available
- ↳ 22 mm advancement with 17 mm notch



Ref. No	Diameter (G)	Length (cm)	Coaxial Needle Length (cm)	With Coaxial Ref. No
GVS1410	14	10	5,5	GVS1410-X
GVS1415	14	15	10,5	GVS1415-X
GVS1420	14	20	15,5	GVS1420-X
GVS1610	16	10	5,5	GVS1610-X
GVS1615	16	15	10,5	GVS1615-X
GVS1620	16	20	15,5	GVS1620-X
GVS1810	18	10	5,5	GVS1810-X
GVS1815	18	15	10,5	GVS1815-X
GVS1820	18	20	15,5	GVS1820-X
GVS1825	18	25		
GVS1830	18	30		
GVS2010	20	10	5,5	GVS2010-X
GVS2015	20	15	10,5	GVS2015-X
GVS2020	20	20	15,5	GVS2020-X

ESTACORE® AUTOMATIC BIOPSY NEEDLE

- ↳ Easy and one hand operated with two stage arming
- ↳ Echogenic and atraumatic needle tip design
- ↳ 22mm advancement with 17mm notch
- ↳ Ergonomic and light weight design
- ↳ Depth marked canula
- ↳ Choice of two firing buttons
- ↳ Coaxial needles are optionally available



Ref. No	Diameter (G)	Length (cm)	Coaxial Needle Length (cm)	With Coaxial Ref. No
GEC1410	14	10	5,5	GEC1410-X
GEC1415	14	15	10,5	GEC1415-X
GEC1420	14	20	15,5	GEC1420-X
GEC1610	16	10	5,5	GEC1610-X
GEC1615	16	15	10,5	GEC1615-X
GEC1620	16	20	15,5	GEC1620-X
GEC1810	18	10	5,5	GEC1810-X
GEC1815	18	15	10,5	GEC1815-X
GEC1820	18	20	15,5	GEC1820-X
GEC1825	18	25		
GEC1830	18	30		
GEC2010	20	10	5,5	GEC2010-X
GEC2015	20	15	10,5	GEC2015-X
GEC2020	20	20	15,5	GEC2020-X



ALFA GUN® REUSABLE BIOPSY DEVICE

- ↳ Light weight and easy grip
- ↳ Two stage cocking by one hand
- ↳ Two optional triggers with high velocity throw
- ↳ Obtained top quality tissue samples quickly and efficiently
- ↳ Compatible with its own disposable needle



Ref. No	Description
GAG01	Alfa Gun Reusable Biopsy Device



ALFA GUN® BIOPSY NEEDLE

Ref. No	Diameter (G)	Length (cm)
GAC1610	16	10
GAC1615	16	15
GAC1620	16	20
GAC1625	16	25
GAC1630	16	30
GAC1810	18	10
GAC1815	18	15
GAC1820	18	20
GAC1825	18	25
GAC1830	18	30

MAXICORE® REUSABLE BIOPSY GUN

- ↳ Trigger is easily activated
- ↳ Preventing accidental firing
- ↳ Easily positioning and removing
- ↳ Sterilized with conventional methods



Ref. No	Description
GBG01	Maxicore Reusable Biopsy Gun



MAXICORE® BIOPSY GUN NEEDLE

- ↳ Echogenic needle tip
- ↳ 17mm sample notch
- ↳ Depth marked canula
- ↳ Color code
- ↳ Compatible with most common reusable guns in the market




Ref. No	Diameter (G)	Length (cm)
GPN1410	14	10
GPN1415	14	15
GPN1420	14	20
GPN1610	16	10
GPN1615	16	15
GPN1620	16	20
GPN1810	18	10
GPN1815	18	15
GPN1820	18	20
GPN1825	18	25
GPN1830	18	30

MAXICORE-M® BIOPSY GUN NEEDLE

- ↳ Echogenic needle tip
- ↳ 17 mm sample notch
- ↳ Depth marked canula
- ↳ Color code
- ↳ Compatible with Bard Magnum® Reusable Gun (*)

(*) Bard Magnum® is a Trade Mark of Bard Inc.




Ref. No	Diameter (G)	Length (cm)
GPN1410-M	14	10
GPN1415-M	14	15
GPN1420-M	14	20
GPN1610-M	16	10
GPN1615-M	16	15
GPN1620-M	16	20
GPN1810-M	18	10
GPN1815-M	18	15
GPN1820-M	18	20
GPN1825-M	18	25
GPN1830-M	18	30

MAXICORE-U® BIOPSY GUN NEEDLE

- ↳ Echogenic needle tip
- ↳ 17 mm sample notch
- ↳ Depth marked canula
- ↳ Color code
- ↳ Compatible with Pro-Mag Ultra® Reusable Gun (*)

(*) Pro-Mag Ultra is a Trade Mark of Angiotech Inc.



Ref. No	Diameter (G)	Length (cm)
GPN1410-U	14	10
GPN1415-U	14	15
GPN1420-U	14	20
GPN1610-U	16	10
GPN1615-U	16	15
GPN1620-U	16	20
GPN1810-U	18	10
GPN1815-U	18	15
GPN1820-U	18	20
GPN1825-U	18	25
GPN1830-U	18	30



SEMI-AUTOMATIC BIOPSY NEEDLE

New
Design

- ↳ Echogenic and atraumatic needle tip design
- ↳ Light weight handle
- ↳ Adjustable sample notch: 10mm-20mm
- ↳ Depth marked canula
- ↳ Color code
- ↳ Co-axial needles are optionally available



Ref. No	Diameter (G)	Length (cm)	Coaxial Needle Length (cm)	With Coaxial Ref. No
GSNA1410	14	10	5	GSNA1410-X
GSNA1415	14	15	10	GSNA1415-X
GSNA1420	14	20	15	GSNA1420-X
GSNA1610	16	10	5	GSNA1610-X
GSNA1615	16	15	10	GSNA1615-X
GSNA1620	16	20	15	GSNA1620-X
GSNA1810	18	10	5	GSNA1810-X
GSNA1815	18	15	10	GSNA1815-X
GSNA1820	18	20	15	GSNA1820-X
GSNA1825	18	25		GSNA1825-X
GSNA2010	20	10	5	GSNA2010-X
GSNA2015	20	15	10	GSNA2015-X
GSNA2020	20	20	15	GSNA2020-X

SEMI-AUTOMATIC BIOPSY NEEDLE

- ↳ Echogenic and atraumatic needle tip design
- ↳ Light weight handle
- ↳ Adjustable sample notch: 10mm-20mm
- ↳ Depth marked canula
- ↳ Color code
- ↳ Co-axial needles are optionally available



Ref. No	Diameter (G)	Length (cm)	Coaxial Needle Length (cm)	With Coaxial Ref. No
GSN1410	14	10	5	GSN1410-X
GSN1415	14	15	10	GSN1415-X
GSN1420	14	20	15	GSN1420-X
GSN1610	16	10	5	GSN1610-X
GSN1615	16	15	10	GSN1615-X
GSN1620	16	20	15	GSN1620-X
GSN1810	18	10	5	GSN1810-X
GSN1815	18	15	10	GSN1815-X
GSN1820	18	20	15	GSN1820-X
GSN1825	18	25		GSN1825-X
GSN2010	20	10	5	GSN2010-X
GSN2015	20	15	10	GSN2015-X
GSN2020	20	20	15	GSN2020-X

CHIBA BIOPSY NEEDLE

- ↳ cm marking
- ↳ Color coded
- ↳ Echogenic tip



Ref. No	Diameter (G)	Length (cm)
GCB1810	18	10
GCB1815	18	15
GCB1820	18	20
GCB2010	20	10
GCB2015	20	15
GCB2020	20	20
GCB2025	20	25
GCB2210	22	10
GCB2215	22	15
GCB2220	22	20
GCB2225	22	25




COAXIAL NEEDLE

- ↳ Light weight and extra thin wall tubing
- ↳ Trocar tip for straight entry
- ↳ Color code for identifying the accurate diameter
- ↳ Accurate depth placement with centimeter marks





Ref. No	Diameter (G)	Length (cm)
GCX1305	13	5
GCX1310	13	10
GCX1315	13	15
GCX1505	15	5
GCX1510	15	10
GCX1515	15	15
GCX1705	17	5
GCX1710	17	10
GCX1715	17	15
GCX1905	19	5
GCX1910	19	10
GCX1915	19	15



BREAST LOCALIZATION NEEDLE

- ↳ Ultra-sharp and atraumatic needle tip design
- ↳ Marking wire with single hook
- ↳ Depth markers



Ref. No	Diameter (G)	Length (cm)
GBL2005	20	5
GBL2007	20	7
GBL2010	20	10
GBL2012	20	12
GBL2015	20	15

PERCUTANEOUS ACCESS NEEDLE

- ↳ Echogenic needle tip
- ↳ Ultra-thin cannula accepts
0.038" guide wire
- ↳ Depth markers
- ↳ Color code



Ref. No	Diameter (G)	Length (cm)
GPA1810	18	10
GPA1815	18	15
GPA1820	18	20
GPA1825	18	25



TUMOR MARKING NEEDLE

w/PEEK® Marker

- Sterile placement needle with preloaded PEEK® marker
- Accurate target visualization in tumors/organs during external radiation therapy
- Excellent visibility of tumor markers under fluoroscopy, X-ray, CT, MRI and ultrasound
- Knurled surface of the marker prevents migration




Ref. No	Needle (G x cm)	PEEK® Marker (mm x mm)
GTP1415	14 x 15	1 x 3, 4, 5
GTP1420	14 x 20	1 x 3, 4, 5
GTP1425	14 x 25	1 x 3, 4, 5
GTP1815	18 x 15	1 x 3, 4, 5
GTP1820	18 x 20	1 x 3, 4, 5
GTP1825	18 x 25	1 x 3, 4, 5
GTP1915	19 x 15	0.8 x 3, 4, 5
GTP1920	19 x 20	0.8 x 3, 4, 5
GTP1925	19 x 25	0.8 x 3, 4, 5

TUMOR MARKING NEEDLE

w/Gold Marker

- Sterile placement needle with preloaded 99.99 % pure Gold marker
- Accurate target visualization in tumors/organs during external radiation therapy
- Excellent visibility of tumor markers under fluoroscopy, X-ray, CT, MRI and ultrasound
- The ergonomic handle design with safety lock provides secure positioning of gold marker
- Knurled surface of the marker prevents migration




Ref. No	Needle (G x cm)	Gold Marker (mm x mm)
GTM1415	14 x 15	1 x 3, 4, 5
GTM1420	14 x 20	1 x 3, 4, 5
GTM1425	14 x 25	1 x 3, 4, 5
GTM1815	18 x 15	1 x 3, 4, 5
GTM1820	18 x 20	1 x 3, 4, 5
GTM1825	18 x 25	1 x 3, 4, 5
GTM1915	19 x 15	0.8 x 3, 4, 5
GTM1920	19 x 20	0.8 x 3, 4, 5
GTM1925	19 x 25	0.8 x 3, 4, 5



BONE MARROW ASPIRATION NEEDLE

- ↳ Atraumatic ultra sharp needle tip design
- ↳ Adjustable needle length
- ↳ Depth marked scala
- ↳ Light weight and ergonomic handle
- ↳ Easy to use with excellent comfort
- ↳ Color code



Ref. No	Diameter (G)	Length (cm)
GBA1450	14	5
GBA1550	15	5
GBA1650	16	5
GBA1850	18	5

BONE MARROW BIOPSY NEEDLE

- ↳ Atraumatic ultra sharp needle tip design
- ↳ Light weight and ergonomic handle
- ↳ Easy to use with excellent comfort



Ref. No	Diameter (G)	Length (cm)
GBB0810	8	10
GBB0815	8	15
GBB0910	9	10
GBB0915	9	15
GBB1110	11	10
GBB1115	11	15
GBB1310	13	10
GBB1315	13	15

DISPOSABLE NEEDLE GUIDES



- ↳ Used in ultrasound guided procedures like as Transrectal prostate biopsies, Transvaginal biopsies and Transvaginal oocyte retrieval
- ↳ The disposable design reduces the risks of cross contamination and infection
- ↳ All guides accept 16G-18G sizes instruments
- ↳ Compatible with well-known endocavity ultrasound probes in the market



****** Please contact with us for the list of needle guides or your custom made request.

TRUS BIOPSY KIT

- Full Set for Transrectal Ultrasound guided prostate biopsy procedures
- Reduced risk of cross-contamination with sterile and single use kit
- Easier and comfortable procedure with all necessary items in one kit
- Offered in a variety of kit configuration
- Custom-made kit configuration is upon your request



Trust Kit Ref. No.	Needle Guide (1 unit)	Reusable Biopsy Gun Needle (1 unit)	Estacore Disposable Biopsy Device (1 unit)	Chiba Needle (1 unit)	Sample Box (12 bottles)	Probe Cover (2 units)
		16Gx20cm 16Gx25cm 16Gx30cm	16Gx20cm 16Gx25cm 16Gx30cm	18Gx20cm 20Gx20cm 22Gx20cm		
		18Gx20cm 18Gx25cm 18Gx30cm	18Gx20cm 18Gx25cm 18Gx30cm			
GTR-A	+	+		+	+	+
GTR-B	+		+	+	+	+
GTR-C	+	+		+		+
GTR-D	+		+	+		+
GTR-E	+				+	+


- * Parts included in kit are shown as "+" on above lines.
- * Please select your kit from the "Ref. No. Column" according to your requirement and notice your needle sizes for biopsy and chiba additionally.
- * Please indicate your model of Needle Guide and USG Transducer when ordering.

OOCYTE ASPIRATION NEEDLE

- ↳ Unique grip design with bevel marker ensures effective control and comfort
- ↳ 17G and 18G Single Lumen Needle maximize the quality of oocytes retrieved along with minimized operating time
- ↳ Echo marked tip gives extreme ultrasound definition
- ↳ Sharp needle tip ensures minimum trauma
- ↳ Silicone tap fits test tubes securely
- ↳ With 3 way stopcock for flushing line is also available




Ref. No	Needle (G x cm)	Aspiration (cm)	Suction Tube (cm)	Flushing Tube w/stopcock
GOA1735	17 x 35	90	40	-
GOA1735-F	17 x 35	90	40	+
GOA1835	18 x 35	90	40	-
GOA1835-F	18 x 35	90	40	+



UROLOGY CATHETERS

LONG DURABLE DJ URETERAL STENT SET

- ↳ Made by Carbothane
- ↳ Suitable for long term implantation
- ↳ Higher kink resistance
- ↳ High radiopacity
- ↳ Depth marker
- ↳ Long pusher option between 35-100 cm
- ↳ Suture on distal tip



Diameter (F)	Length (cm)	Ref. No: (One End Open)	Ref. No: (Both End Open)
3	12	GDJ312-OEC-CAR	GDJ312-BEO-CAR
	14	GDJ314-OEC-CAR	GDJ314-BEO-CAR
	16	GDJ316-OEC-CAR	GDJ316-BEO-CAR
	20	GDJ320-OEC-CAR	GDJ320-BEO-CAR
4	12	GDJ412-OEC-CAR	GDJ412-BEO-CAR
	14	GDJ414-OEC-CAR	GDJ414-BEO-CAR
	16	GDJ416-OEC-CAR	GDJ416-BEO-CAR
	20	GDJ420-OEC-CAR	GDJ420-BEO-CAR
	22	GDJ422-OEC-CAR	GDJ422-BEO-CAR
	24	GDJ424-OEC-CAR	GDJ424-BEO-CAR
4.7	26	GDJ426-OEC-CAR	GDJ426-BEO-CAR
	12	GDJ4712-OEC-CAR	GDJ4712-BEO-CAR
	16	GDJ4716-OEC-CAR	GDJ4716-BEO-CAR
	20	GDJ4720-OEC-CAR	GDJ4720-BEO-CAR
	22	GDJ4722-OEC-CAR	GDJ4722-BEO-CAR
	24	GDJ4724-OEC-CAR	GDJ4724-BEO-CAR
4.8	26	GDJ4726-OEC-CAR	GDJ4726-BEO-CAR
	28	GDJ4728-OEC-CAR	GDJ4728-BEO-CAR
	12	GDJ4812-OEC-CAR	GDJ4812-BEO-CAR
	20	GDJ4820-OEC-CAR	GDJ4820-BEO-CAR
	22	GDJ4822-OEC-CAR	GDJ4822-BEO-CAR
	24	GDJ4824-OEC-CAR	GDJ4824-BEO-CAR
5	26	GDJ4826-OEC-CAR	GDJ4826-BEO-CAR
	28	GDJ4828-OEC-CAR	GDJ4828-BEO-CAR
	16	GDJ516-OEC-CAR	GDJ516-BEO-CAR
	22	GDJ522-OEC-CAR	GDJ522-BEO-CAR
	24	GDJ524-OEC-CAR	GDJ524-BEO-CAR
6	26	GDJ526-OEC-CAR	GDJ526-BEO-CAR
	28	GDJ528-OEC-CAR	GDJ528-BEO-CAR
	20	GDJ620-OEC-CAR	GDJ620-BEO-CAR
	22	GDJ622-OEC-CAR	GDJ622-BEO-CAR
	24	GDJ624-OEC-CAR	GDJ624-BEO-CAR
	26	GDJ626-OEC-CAR	GDJ626-BEO-CAR
	28	GDJ628-OEC-CAR	GDJ628-BEO-CAR
	30	GDJ630-OEC-CAR	GDJ630-BEO-CAR
8	24	GDJ724-OEC-CAR	GDJ724-BEO-CAR
	26	GDJ726-OEC-CAR	GDJ726-BEO-CAR
	28	GDJ728-OEC-CAR	GDJ728-BEO-CAR
	30	GDJ730-OEC-CAR	GDJ730-BEO-CAR
	22	GDJ822-OEC-CAR	GDJ822-BEO-CAR
	24	GDJ824-OEC-CAR	GDJ824-BEO-CAR
	26	GDJ826-OEC-CAR	GDJ826-BEO-CAR
	28	GDJ828-OEC-CAR	GDJ828-BEO-CAR
	30	GDJ830-OEC-CAR	GDJ830-BEO-CAR

↳ Set includes DJ stent + Pusher + Guide wire + Two Clamps

DOUBLE J URETERAL STENT SET

- ↳ Soft polyurethane material
- ↳ Higher kink resistance
- ↳ High radiopacity
- ↳ Depth marker
- ↳ Long pusher option between 35 cm - 100 cm
- ↳ Suture on distal tip




DJ + PUSHER ONLY

Please Use Ref. No as GDJPXX-OEC or BE0

HYDROPHILIC DJ STENT SET

Please Use Ref. No as GDJXX-HOEC or HBE0



Diameter (F)	Length (cm)	Ref. No: (One End Open)	Ref. No: (Both End Open)
3	12	GDJ312-OEC	GDJ312-BE0
	14	GDJ314-OEC	GDJ314-BE0
	16	GDJ316-OEC	GDJ316-BE0
	20	GDJ320-OEC	GDJ320-BE0
4	12	GDJ412-OEC	GDJ412-BE0
	14	GDJ414-OEC	GDJ414-BE0
	16	GDJ416-OEC	GDJ416-BE0
	20	GDJ420-OEC	GDJ420-BE0
	22	GDJ422-OEC	GDJ422-BE0
	24	GDJ424-OEC	GDJ424-BE0
4.7	26	GDJ426-OEC	GDJ426-BE0
	12	GDJ4712-OEC	GDJ4712-BE0
	16	GDJ4716-OEC	GDJ4716-BE0
	20	GDJ4720-OEC	GDJ4720-BE0
	22	GDJ4722-OEC	GDJ4722-BE0
	24	GDJ4724-OEC	GDJ4724-BE0
Adjustable	26	GDJ4726-OEC	GDJ4726-BE0
	28	GDJ4728-OEC	GDJ4728-BE0
	22/30	VDJ47-OEC	VDJ47-BE0
	12	GDJ4812-OEC	GDJ4812-BE0
4.8	16	GDJ4816-OEC	GDJ4816-BE0
	20	GDJ4820-OEC	GDJ4820-BE0
	22	GDJ4822-OEC	GDJ4822-BE0
	24	GDJ4824-OEC	GDJ4824-BE0
	26	GDJ4826-OEC	GDJ4826-BE0
	28	GDJ4828-OEC	GDJ4828-BE0
Adjustable	22/30	VDJ48-OEC	VDJ48-BE0
5	16	GDJ516-OEC	GDJ516-BE0
	20	GDJ520-OEC	GDJ520-BE0
	22	GDJ522-OEC	GDJ522-BE0
	24	GDJ524-OEC	GDJ524-BE0
	26	GDJ526-OEC	GDJ526-BE0
	28	GDJ528-OEC	GDJ528-BE0
6	20	GDJ620-OEC	GDJ620-BE0
	22	GDJ622-OEC	GDJ622-BE0
	24	GDJ624-OEC	GDJ624-BE0
	26	GDJ626-OEC	GDJ626-BE0
	28	GDJ628-OEC	GDJ628-BE0
	30	GDJ630-OEC	GDJ630-BE0
Adjustable	22/30	VDJ06-OEC	VDJ06-BE0
7	22	GDJ722-OEC	GDJ722-BE0
	24	GDJ724-OEC	GDJ724-BE0
	26	GDJ726-OEC	GDJ726-BE0
	28	GDJ728-OEC	GDJ728-BE0
	30	GDJ730-OEC	GDJ730-BE0
	22	GDJ822-OEC	GDJ822-BE0
8	24	GDJ824-OEC	GDJ824-BE0
	26	GDJ826-OEC	GDJ826-BE0
	28	GDJ828-OEC	GDJ828-BE0
	30	GDJ830-OEC	GDJ830-BE0

↳ Set includes DJ stent + Pusher + Guide wire + Two Clamp

URETERAL CATHETER

- ↳ Wide size range
- ↳ Depth marker
- ↳ Luer lock adapter
- ↳ Metal Mandren

URETERAL CATHETER CLOSED TIP		
Ref. No	Diameter (F)	Length (cm)
GUC370-CE	3	70
GUC470-CE	4	70
GUC570-CE	5	70
GUC670-CE	6	70
GUC770-CE	7	70

URETERAL CATHETER OPEN TIP		
Ref. No	Diameter (F)	Length (cm)
GUC370-OE	3	70
GUC470-OE	4	70
GUC570-OE	5	70
GUC670-OE	6	70
GUC770-OE	7	70

SINGLE J URETERAL STENT SET

- ↳ Two different colors single J catheters for left and right side
- ↳ Guide wire
- ↳ Dilators 7 F to 9 F
- ↳ Urine Bag Connector



Ref. No	Diameter (F)	Length (cm)
GSJ 670	6	70
GSJ 770	7	70
GSJ 870	8	70

RE-ENTRY MALECOT CATHETER

- ↳ Two flexible stylets
- ↳ High radiopacity
- ↳ Smooth surface
- ↳ Drainage bag adaptor

Ref. No	Diameter (F)	Body Length (cm)	Tail Length (cm)
GRM 1430	14	30	15
GRM 1630	16	30	15
GRM 1830	18	30	15
GRM 2030	20	30	15

PERCUTANEOUS NEPHROSTOMY SET



- ↳ Pigtail or malecot tip catheter with:
- ↳ Puncture needle
- ↳ Dilators
- ↳ Stiff guide wire or lunderquist guide wire options
- ↳ Scalpel
- ↳ Drainage bag connector

NEPHROSTOMY SET

Ref. No	Diameter (F)	With Locking	Hydrophilic	Hydrophilic & Locking
GPNS06	6	GPNS06-L	GPNSH06	GPNSH06-L
GPNS08	8	GPNS08-L	GPNSH08	GPNSH08-L
GPNS10	10	GPNS10-L	GPNSH10	GPNSH10-L
GPNS12	12	GPNS12-L	GPNSH12	GPNSH12-L
GPNS14	14	GPNS14-L	GPNSH14	GPNSH14-L
GPNS16	16	GPNS16-L	GPNSH16	GPNSH16-L

NEPHROSTOMY CATHETER

Ref. No	Diameter (F)	With Locking	Hydrophilic	Hydrophilic & Locking
GNC06	6	GNC06-L	GNCH06	GNCH06-L
GNC08	8	GNC08-L	GNCH08	GNCH08-L
GNC10	10	GNC10-L	GNCH10	GNCH10-L
GNC12	12	GNC12-L	GNCH12	GNCH12-L
GNC14	14	GNC14-L	GNCH14	GNCH14-L
GNC16	16	GNC16-L	GNCH16	GNCH16-L

DIRECT PUNCTURE NEPHROSTOMY SET

- ↳ Pigtail nephrostomy catheter with metal stiffener and trocar point stylet
- ↳ Easy penetration with sharp beveled trocar and tapered pigtail catheter tip
- ↳ Large drainage holes
- ↳ Hydrophilic coated and locked catheter types are optional



Ref. No	Description	Size (Fr)	Lenght (cm)
GDPN08	Direct Puncture Nephrostomy Set	8	30
GDPN10	Direct Puncture Nephrostomy Set	10	30
GDPN12	Direct Puncture Nephrostomy Set	12	30
GDPN14	Direct Puncture Nephrostomy Set	14	30
GDPN08-L	Direct Puncture Nephrostomy Set with Locking	8	30
GDPN10-L	Direct Puncture Nephrostomy Set with Locking	10	30
GDPN12-L	Direct Puncture Nephrostomy Set with Locking	12	30
GDPN14-L	Direct Puncture Nephrostomy Set with Locking	14	30
GDPNH08	Direct Puncture Nephrostomy Set Hydrophilic Coated	8	30
GDPNH10	Direct Puncture Nephrostomy Set Hydrophilic Coated	10	30
GDPNH12	Direct Puncture Nephrostomy Set Hydrophilic Coated	12	30
GDPNH14	Direct Puncture Nephrostomy Set Hydrophilic Coated	14	30
GDPNH08-L	Direct Puncture Nephrostomy Set with Locking and Hydrophilic Coted	8	30
GDPNH10-L	Direct Puncture Nephrostomy Set with Locking and Hydrophilic Coted	10	30
GDPNH12-L	Direct Puncture Nephrostomy Set with Locking and Hydrophilic Coted	12	30
GDPNH14-L	Direct Puncture Nephrostomy Set with Locking and Hydrophilic Coted	14	30

URETERAL BALLOON DILATOR




- ↳ Used for dilatation of ureter
- ↳ 2 ways catheter
- ↳ Platinum markers at both end
- ↳ 18 to 20 atm rated burst pressure
- ↳ Supplied with stopcock

Ref. No	Balloon Diameter (mm)	Balloon Length (cm)	Shaft (F)	Usable Length (cm)	Recommended Guide wire (inch)	Balloon RBP (atm)
GBD0544	4	4	5	70	.38"	20
GBD0546	4	6	5	70	.38"	20
GBD5410	4	10	5	70	.38"	20
GBD0554	5	4	5	70	.38"	20
GBD0556	5	6	5	70	.38"	20
GBD5510	5	10	5	70	.38"	18
GBD0564	6	4	5	70	.38"	18
GBD0566	6	6	5	70	.38"	18
GBD5610	6	10	5	70	.38"	18
GBD0444	4	4	4	70	.18"	20
GBD0454	5	4	4	70	.18"	20

SUPRAPUBIC DRAINAGE CATHETER SET

- ↳ Silicone coated Foley Catheter
- ↳ Trocar with peel away sheath
- ↳ Scalpel



Ref. No	Catheter Diameter (F)	Catheter Length (cm)
GSP10	10	35
GSP12	12	35
GSP14	14	35
GSP16	16	35
GSP18	18	35



AMPLATZ RENAL DILATOR SET

- ↳ Used for renal dilatation before stone removal procedures
- ↳ High radioopacity
- ↳ Wide range of dilator options: 8F to 30F
- ↳ 4 Pcs sheath and 70 cm introducer catheter



Ref. No	Description	Size
GUDS-12	Amplatz Renal Dilator Set	3F - 30 F

STONE BASKET

- ↳ Stainless steel and nitinol alternatives
- ↳ Designed for endoscopic manipulation
- ↳ Detachable handle



STAINLESS STEEL

Ref. No	Diameter (F)	Length (cm)	BASKET TYPE	WIRE TYPE
GSBH390	3	90	HELICAL	Round
GSBS390			PARALLEL	Flat
GSBH3110	3	110	HELICAL	Round
GSBS3110			PARALLEL	Flat
GSBH490	4	90	HELICAL	Round
GSBS490			PARALLEL	Flat
GSBH4110	4	110	HELICAL	Round
GSBS4110			PARALLEL	Flat
GSBH590	5	90	HELICAL	Round
GSBS590			PARALLEL	Flat
GSBH5110	5	110	HELICAL	Round
GSBS5110			PARALLEL	Flat

NITINOL


Ref. No	Diameter (G)	Length (cm)	BASKET TYPE	WIRE TYPE	Zero Tip Ref. No
GSBN390	1.9	90	HELICAL	Round	GSBZ1990
GSBN3110	1.9	110	HELICAL	Round	GSBZ19110
GSBN3110	1.9	120	HELICAL	Round	GSBZ19120
GSBN490	2.2	90	HELICAL	Round	GSBZ2290
GSBN4110	2.2	110	HELICAL	Round	GSBZ22110
GSBN390	3	90	HELICAL	Round	GSBZ390
GSBN3110	3	110	HELICAL	Round	GSBZ3110
GSBN490	4	90	HELICAL	Round	GSBZ490
GSBN4110	4	110	HELICAL	Round	GSBZ4110
GSBN590	5	90	HELICAL	Round	GSBZ590
GSBN5110	5	110	HELICAL	Round	GSBZ5110

STONE GRASPER

- ↳ Stainless steel
- ↳ Designed for endoscopic manipulation
- ↳ Detachable handle



Ref. No	Description	Diameter (F)	Length (cm)
GSG390	Grasping Forceps, Three Pong	3	90
GSG490	Grasping Forceps, Three Pong	4	90
GSG590	Grasping Forceps, Three Pong	5	90



NEPHROSTOMY BALLOON DILATOR

- ↳ Used for dilatation of renal tract
- ↳ 2 ways catheter
- ↳ Platinum markers at both end
- ↳ 17 atm rated burst pressure
- ↳ Supplied with 30 F PTFE Sheath



Ref. No	Balloon Diameter (mm)	Balloon Length (cm)	Shaft (F)	Usable Length (cm)	Recommended Guide wire (inch)	Balloon RBP (atm)
GNBD10	10	12	7	55	.38"	17



URETERAL ACCESS SHEATH

↳ Ureteral access sheath provides ureteral dilation and continuous working channel for the introduction of endoscopes and instruments during ureteral access procedures

Ref. No	Diameter (F)	Length (cm)
GUA1125	9.5/11.5	25
GUA1138	9.5/11.5	38
GUA1145	9.5/11.5	45
GUA1225	10/12	25
GUA1238	10/12	38
GUA1245	10/12	45
GUA1425	12/14	25
GUA1438	12/14	38
GUA1445	12/14	45

NOTES



GEOTEK®

Healthcare Products

İvedik OSB Ağaç-Metal San. Sit.

1436. Sokak No: 12

06370 Ostim - Ankara / TÜRKİYE

Tel: +90 (312) 395 98 52

Fax: +90 (312) 395 98 53

e-mail: info@geotekmedical.com

www.geotekmedical.com